

# SensiFAST™ SYBR & Fluorescein One-Step Kit

Shipping: On Dry/Blue Ice Catalog Numbers

Exp. Date: See vial BIO-75001: 100 x 20µl reactions: 1 x 1ml

Batch No.: See vial BIO-75005: 500 x 20µl reactions: 5 x 1ml

Concentration: see vial



A Meridian Life Science® Company

Store at -20°C

## Storage and Stability:

The SensiFAST SYBR & Fluorescein One-Step Kit is shipped on Dry/Blue Ice. All kit components should be stored at -20°C upon receipt. Excessive freeze/thawing is not recommended. When stored under optimum conditions, the reagents are stable for a minimum of 12 months from date of purchase.

## Quality Control:

Bioline operates under ISO 9001 Management System. The SensiFAST SYBR & Fluorescein One-Step Kit and its components are extensively tested for activity, processivity, efficiency, heat activation, sensitivity and absence of nuclease contamination and nucleic acid contamination.

## Safety Precautions:

Harmful if swallowed. Irritating to eyes, respiratory system and skin. Please refer to the material safety data sheet for further information.

## Notes:

For research use only.

## Description

The SensiFAST™ SYBR & Fluorescein One-Step Kit has been formulated for highly reproducible first-strand cDNA synthesis and subsequent real-time PCR in a single tube. A combination of the latest advances in buffer chemistry together with a reverse transcriptase and hot-start DNA polymerase system ensures that SensiFAST SYBR & Fluorescein One-Step Kit produces fast, highly-specific and ultra-sensitive one-step RT-qPCR.

The SensiFAST SYBR & Fluorescein One-Step Kit consists of a 2x SensiFAST SYBR One-Step mix, as well as separate reverse transcriptase and RiboSafe RNase Inhibitor.

## Kit components

| Reagent   | 100 x 20µl Reactions | 500 x 20µl Reactions |
|---|----------------------|----------------------|
| SensiFAST™ SYBR & Fluorescein One-Step mix (2x) | 1 x 1ml              | 5 x 1ml              |
| RiboSafe RNase Inhibitor                        | 1 x 40µl             | 1 x 200µl            |
| Reverse transcriptase                           | 1 x 20µl             | 1 x 100µl            |
| DEPC-H <sub>2</sub> O                           | 1 x 1.8ml            | 2 x 1.8ml            |

## Instrument compatibility

SensiFAST SYBR & Fluorescein One-Step Kit has been optimized for use in SYBR Green-based real-time RT-PCR on the real-time instruments listed in the following compatibility table, each of these instruments having the capacity to analyze the real-time PCR data with the passive reference signal either on or off. The kit is also compatible with several instruments that do not require the use of fluorescein, such as the Qiagen (Corbett) Rotor-Gene™ 6000, the Bio-Rad CFX96 or the Roche LightCycler® 480.

| Manufacturer | Model                 |
|--------------|-----------------------|
| Bio-Rad      | iCycler™, MyiQ™, iQ™5 |

## General considerations

When handling RNA, it is important to use RNase-free plasticware and reagents. We also recommend performing RNA work in an isolated area. To help prevent any carry-over DNA contamination we recommend that separate areas are maintained for reaction set-up, PCR amplification and any post-PCR gel analysis. It is essential that any tubes containing amplified PCR product are not opened in the PCR set-up area.

**Primers:** The sequence and concentration of the primers, as well as amplicon length, can be critical for specific amplification, yield and overall efficiency of any RT-qPCR. We strongly recommend taking the following points into consideration when designing and running your RT-qPCR:

- use primer-design software, such as Primer3 or visual OMP™ (<http://frodo.wi.mit.edu/primer3/> and DNA Software, Inc. <http://dnasoftware.com/>, respectively). Primers should have a melting temperature (T<sub>m</sub>) of approximately 60°C
- optimal amplicon length should be 80-200bp, and should not exceed 400bp
- final primer concentration of 400nM is suitable for most SYBR-Green based reactions, however to determine the optimal concentration we recommend titrating in the range 0.1-1µM
- use an equimolar primer concentration
- where possible, use intron-spanning primers to avoid amplification from genomic DNA

**Template:** It is important that the RNA template is intact and devoid of DNA or contaminating inhibitors of both reverse transcription and PCR. For high purity RNA, we recommend using the Bioline ISOLATE RNA Mini Kit (BIO-52043). RNA stocks and dilutions should be made in DEPC-treated Water (BIO-38030) to avoid any RNase-mediated degradation.

The recommended amount of template for one-step RT-qPCR is dependent upon the type of RNA used.

- **total RNA:** purified total RNA can be used in the range from 1pg to 1µg per 20µl reaction
- **mRNA:** purified mRNA can be used from 0.01pg per 20µl reaction

**MgCl<sub>2</sub>:** The MgCl<sub>2</sub> concentration in the 1x reaction mix is 3mM. In the majority of RT-qPCR conditions this is optimal for both the reverse transcriptase and the hot-start DNA polymerase. If necessary, we suggest titrating the MgCl<sub>2</sub> to a maximum of 5mM.

**RT-PCR controls:** It is important to detect the presence of contaminating DNA that may affect the reliability of the data. Always include a no-RT control by omitting the reverse transcriptase from the reaction.

**Optional Fluorescein well-factor correction:** SYBR Fluorescein Kit is premixed with fluorescein, so that fluorescence emitted by fluorescein can be optionally detected on certain real-time instruments. If your real-time instrument has the capability of using fluorescein and you wish to use this option, then this option must be selected by the user in the software (see *notice to purchaser No. 5 in Trademark and Licensing Information*).

## Procedure

**Reaction mix composition:** Prepare a RT-PCR mastermix. The volumes given below are based on a standard 20µl final reaction mix and can be scaled accordingly.

| Reagent                                       | Volume     | Final concentration |
|---|------------|---------------------|
| 2x SensiFAST™ SYBR & fluorescein One-Step Mix | Flu- 10µl  | 1x                  |
| 10µM Forward Primer                           | 0.8µl      | 400nM               |
| 10µM Reverse Primer                           | 0.8µl      | 400nM               |
| Reverse transcriptase                         | 0.2µl      | -                   |
| RiboSafe RNase Inhibitor                      | 0.4µl      | -                   |
| H <sub>2</sub> O                              | up to 16µl |                     |
| Template                                      | 4µl        |                     |
| <b>20µl Final volume</b>                      |            |                     |

**Suggested RT-qPCR conditions:** The following RT-qPCR conditions are suitable for the SensiFAST SYBR & Fluorescein One-Step Kit with the majority of amplicons and real-time PCR instruments. However, the cycling conditions can be varied to suit different machine-specific protocols. SensiFAST SYBR & Fluorescein One-Step Kit is compatible with either 3-step or 2-step cycling:

### • 3-step cycling

| Cycles | Temperature          | Time            | Notes   |
|--------|----------------------|-----------------|---|
| 1      | 45°C                 | 10min           | Reverse transcription   |
| 1      | 95°C                 | 2min            | Polymerase activation   |
| 40     | 95°C<br>60°C<br>72°C | 5s<br>10s<br>5s | Denaturation<br>Annealing<br>Extension (acquire at end of step) |

### • 2-step cycling

| Cycles | Temperature  | Time      | Notes  |
|--------|--------------|-----------|--|
| 1      | 45°C         | 10min     | Reverse transcription  |
| 1      | 95°C         | 2min      | Polymerase activation  |
| 40     | 95°C<br>60°C | 5s<br>20s | Denaturation<br>Annealing/extension (acquire at end of step) |

**Optional analysis:** After the reaction has reached completion, refer to the instrument instructions for the option of melt-profile analysis.

## Troubleshooting Guide

| Problem  | Possible Cause   | Recommendation  |
|--|--|---|
| No amplification trace<br>AND<br>No product on agarose gel | Activation time too short  | Ensure SensiFAST SYBR One-Step mix is activated for a minimum of 2min at 95°C before cycling  |
|  | Error in protocol setup  | Verify that correct reagent concentrations, volumes, dilutions and storage conditions have been used  |
|  | Suboptimal primer design   | Use primer design software or validated assay. Test assay on a control template   |
|  | Incorrect concentration of primers   | Use primer concentrations between 100nM and 1µM   |
|  | Template degraded  | Re-isolate your template from the sample material or use freshly prepared template dilution. We recommend using the ISOLATE RNA kits for template preparation and DEPC-treated water for resuspension or dilution of the template |
|  | Primers degraded   | Use newly synthesized primers   |
|  | Template contaminated with RT-PCR inhibitors   | Further dilute template before RT-PCR or purify template and resuspend it in DEPC-treated water   |
|  | Template concentration too low   | Increase concentration used   |
| Cycling conditions not optimal                             | Increase extension/annealing time, increase cycle number, reduce annealing temperature |   |

## Troubleshooting Guide (Continued)

| Problem   | Possible Cause  | Recommendation   |
|---|---|--|
| No amplification trace<br>AND<br>PCR product present on agarose gel | Error in instrument setup                                   | Check that the acquisition settings are correct during cycling   |
| Non-specific amplification product<br>AND<br>Primer-dimers          | Inefficient reverse transcription                           | Extend reverse transcription time up to 20min and/or increase the temperature up to 48°C   |
|   | Suboptimal primer design                                    | Redesign primers using appropriate software, or use validated primers  |
|   | Primer concentration too high                               | Test dilution series of primer concentrations until primer-dimer/non-specific amplification products disappear   |
|   | Primer concentration too low                                | Use primer concentration between 100nM and 1µM   |
|   | Primer annealing temperature too low                        | Increase PCR annealing temperature up to 65°C or until primer-dimer/non-specific amplification products disappear                                      |
|   | Template concentration too low                              | Increase template concentration  |
|   | Template concentration too high                             | Reduce template concentration until non-specific products disappear  |
| Variability between replicates                                      | Extension time too long                                     | Reduce extension time to determine whether non-specific products are reduced   |
|   | Error in reaction set-up                                    | Prepare large volume mastermix   |
| Late amplification trace  | Air bubbles in reaction mix                                 | Centrifuge reaction samples/plate prior to running on a real-time instrument   |
|   | Inefficient reverse transcription                           | Extend reverse transcription time up to 20min and/or increase the temperature up to 48°C   |
|   | Activation time too short                                   | Ensure SensiFAST SYBR & Fluorescein One-Step mix is activated for a minimum of 1min at 95°C before cycling   |
|   | Annealing temperature too high                              | Decrease annealing temperature in steps of 2°C   |
|   | Extension time too short                                    | Double extension time to determine whether the cycle threshold (C <sub>T</sub> ) is affected   |
|   | Template concentration too low                              | Increase concentration if possible   |
|   | Template is degraded  | Re-isolate template from sample material or use freshly prepared template dilution   |
|   | Suboptimal primer design                                    | Redesign primers using appropriate software, or use validated primers  |
|   | Primer concentration too low                                | Increase concentration of primers in 100nM increments  |
| RNase contamination   | Ensure RNase inhibitor is added before addition of template |  |
| PCR efficiency below 90%  | Extension time too short                                    | Increase extension time  |
|   | Primer concentration too low                                | Increase concentration of primers in 100nM increments  |
|   | Suboptimal primer design                                    | Redesign primers using appropriate software or use validated primers   |
| PCR efficiency above 110%   | Template is degraded or contains PCR inhibitors             | Re-isolate template from sample material, or use freshly prepared template dilution, or purify template and resuspend it in water                      |
|   | Non-specific amplification and/or primer-dimers             | Use 4% agarose gel electrophoresis to confirm presence of non-specific amplification products. See above for preventing/removing non-specific products |

## Associated Products

| Product                      | Description  | Pack Size                         | Cat No.                             |
|------------------------------|--|-----------------------------------|-------------------------------------|
| ISOLATE Genomic DNA Mini kit | Rapid isolation of DNA from a variety of samples   | 10 Preps<br>50 Preps<br>250 Preps | BIO-52031<br>BIO-52032<br>BIO-52033 |
| ISOLATE Plant DNA Mini kit   | Rapid isolation of DNA from a variety of plant samples   | 10 Preps<br>50 Preps<br>250 Preps | BIO-52034<br>BIO-52035<br>BIO-52036 |
| ISOLATE RNA Mini Kit         | Fast and efficient isolation of extremely pure total RNA from a variety of samples                 | 10 Preps<br>50 Preps<br>250 Preps | BIO-52039<br>BIO-52040<br>BIO-52041 |
| ISOLATE Plant RNA Mini Kit   | Fast and efficient isolation of extremely pure total RNA from a variety of plant samples           | 10 Preps<br>50 Preps<br>250 Preps | BIO-52042<br>BIO-52043<br>BIO-52044 |
| TRIsure™                     | Quick isolation of high-quality RNA from a variety of sources for subsequent use in cDNA synthesis | 100ml<br>200ml                    | BIO-38032<br>BIO-38033              |
| cDNA Synthesis Kit           | Fully optimized to generate maximum yields of full-length cDNA from RNA                            | 30 Reactions<br>100 Reactions     | BIO-65025<br>BIO-65026              |
| Agarose                      | Molecular biology grade agarose  | 100g<br>500g                      | BIO-41026<br>BIO-41025              |
| PCR Water                    | Ultra-pure (18.2MΩ) molecular biology grade water  | 10 x 10ml                         | BIO-37080                           |
| DEPC-treated Water           | Deionized, high-quality molecular grade water treated with DEPC. Ideal for use in all RNA work     | 10 x 10ml<br>1 Liter              | BIO-38030<br>BIO-38031              |

## Technical Support

If the troubleshooting guide does not solve the difficulty you are experiencing, please contact Technical Support with details of reaction setup, cycling conditions and relevant data.

Email: [tech@bioline.com](mailto:tech@bioline.com)

## Trademark and Licensing Information

1). Trademarks: SensiFAST™ (Bioline Reagents Ltd), SYBR® (Molecular Probes), ROX™, LightCycler™ (Roche), StepOne™ (ABI), RotorGene™ (Qiagen), LightCycler® (Roche), iCycler™ MyiQ™, IQ™ (Bio-Rad).

2). Purchase of this product includes limited right to use the supplied amount of SYBR® Green I Stain patented by Molecular Probes, Inc.

3) Notice to Purchaser: Limited License. Use of this product may be covered by one or more of the following US patents: 6,127,155, 5,677,152 (claims 1 to 23 only), 5,773,258 (claims 1 and 6 only). The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

4) SensiFAST products are manufactured by Bioline Reagents Ltd.

5) Notice to Purchaser: No rights are conveyed with respect to US patent 5,928,907

6) Notice to Purchaser: Licensed under US patents 5,338,671 and 5,587,287 and corresponding patents in other countries

Bioline Reagents Ltd  
UNITED KINGDOM  
  
Tel: +44(0)20 8830 5300  
Fax: +44 (0)20 8452 2822

Bioline USA Inc.  
USA  
  
Tel: +1 508 880 8990  
Fax: +1 508 880 8993

Bioline GmbH  
GERMANY  
  
Tel: +49(0)33 7168 1229  
Fax: +49 (0)337168 1244

Bioline (Aust) Pty. Ltd  
AUSTRALIA  
  
Tel: +61 (0)2 9209 4180  
Fax: +61 (0)2 9209 4763