


SPHERO™ AccuCount Particles

- Contains particles with a known number of particles per mL
- Easy to use and cost effective
- Available in several sizes to accommodate the cell size to be counted
- Provided as blank, fluorescent, Rainbow, or Ultra Rainbow particles that can be used in multiple fluorescent channels of flow cytometers

The SPHERO™ AccuCount Particles are designed to be used as reference particles with known number of particles per mL for counting the absolute cell number by flow cytometry. The SPHERO™ AccuCount Particles are very easy to use and are cost effective. The AccuCount Fluorescent Particles are fluorescent in FITC, PE and PE-Cy5 channel. Both AccuCount Fluorescent and AccuCount Blank (nonfluorescent) Particles are available in various particle sizes to accommodate the size of the cells to be counted. In addition, Spherotech also manufactures the AccuCount Rainbow Fluorescent Particles and AccuCount Ultra Rainbow Fluorescent Particles for detection in more fluorescent channels.

Figure 48 Histograms of Cat. No. ACFP-100-3 (AccuCount Fluorescent Particles, 10⁶/mL, 10.2 μm, 3 mL) on a Stratadigm S1400



Particle Type and Surface	Size, μm	Catalog No.	Unit
AccuCount Blank, 10 ⁶ /mL	2.0-2.4	ACBP-20-10	10 mL
AccuCount Blank, 10 ⁶ /mL	5.0-5.9	ACBP-50-10	10 mL
AccuCount Blank, 10 ⁶ /mL	7.0-7.9	ACBP-70-10	10 mL
AccuCount Blank, 10 ⁶ /mL	8.0-12.9	ACBP-100-10	10 mL
AccuCount Fluorescent, 10 ⁶ /mL	5.0-5.9	ACFP-50-5	5 mL
AccuCount Fluorescent, 10 ⁶ /mL	7.0-7.9	ACFP-70-5	5 mL
AccuCount Fluorescent, 10 ⁶ /mL	7.0-7.9	ACFP-70-10	10 mL
AccuCount Fluorescent, 10 ⁶ /mL	8.0-12.9	ACFP-100-3	3 mL
AccuCount Rainbow Fluorescent, 10 ⁶ /mL	8.0-12.9	ACRFP-100-3	3 mL
AccuCount Ultra Rainbow Fluorescent, 10 ⁶ /mL	5.0-5.9	ACURFP-50-10	10 mL
AccuCount Fluorescent Particle Kit. Contains 4 bottles: 0.5x10 ⁶ /mL, 1.0x10 ⁶ /mL, 0.5x10 ⁷ /mL, & 1.0x10 ⁷ /mL	5.0-5.9	ACFP-50-4K	4x1 mL

Figure 49 Histograms of 5 micron fluorescent particles spiked with Cat. No. ACFP-100-3




Figure 50 Histograms of Cat. No. ACBP-100-3 (AccuCount Blank Particles, 10⁶/mL, 10.2 μm, 10 mL)

